

[illegible]

First Wash Protocol: (wash-1.chm) To wash any trace of Piperidine from resin.

1. Dispense System Fluid DMF [1] 4000 µl to RVBlock[1-2]
2. Mix RVBlock for 3 minutes at 450 rpm
3. Empty RVBlock for 2 minutes
- 4.
5. Repeat from step 1, 8 times
6. Return

Couple Protocol: (couple.chm) HOBt/DIC double couple (Do not use repeats or returns in this file)

1. Dispense System Fluid DMF [1] 300 µl to RVBlock[1-2]
(the porous frit will soak up this dead volume)
2. <Dispense Matrix> (Dispense List)
—this is a place holder in software for dispensing amino acids; amino acids are dissolved in 0.5 M HOBt
3. Transfer **1200** µl from DIC [1][DIC] to RVBlock[1-2] using DMF [1]
** Advanced—recommend always using a 250 µl overdraw for Transfer steps, other parameters as default. Wash ratio=2
4. Mix RVBlock for 1 hour at 450 rpm
5. Empty RVBlock for 2 minutes
- 6.
7. Dispense System Fluid DMF [1] 300 µl to RVBlock[1-2]
8. <Dispense Matrix>
9. Transfer **1200** µl from DIC [1][DIC] to RVBlock[1-2] using DMF [1]
** Advanced—recommend always using a 250 µl overdraw for Transfer steps. Wash ratio=2
10. Mix RVBlock for 1 hour at 450 rpm
11. Empty RVBlock for 2 minutes

Note: **1200** is the amount from above **Calculation**, each synthesis will be different.

Second Wash Protocol: (wash-2.chm) To wash any remaining amino acid that was not coupled.

1. Dispense System Fluid DMF [1] 4000 µl to RVBlock[1-2]
2. Mix RVBlock for 3 minutes at 450 rpm
3. Empty RVBlock for 2 minutes
- 4.
5. Repeat from step 1, 3 times
6. Return

Final-Deprotection Protocol: (fin-depr.chm)

1. Goto Chemfile C:\ACT\DEP-2.CHM, line 1
2. Goto Chemfile C:\ACT\WASH-1.CHM, line 1
- 3.
4. Return

DCM-FL.CHM:

1. Flush Arm 1 with DCM [1]
2. Repeat from step 1, 6 times
- 3.
4. Dispense System Fluid DCM [1] 4000µl to RVBlock[1-2]
5. Mix for 3.00 minutes at 450rpm
6. Empty RVBlock for 2 minutes
- 7.
8. Repeat from step 4, 8 times
9. Empty RVBlock for 20 minutes

10.

CLEAVAGE.CHM:

1. Empty CleavageBlock for 2.0 minutes
2. Transfer 4000µl from CocktailSolution[1][] to CleavageBlock[1-2] using DCM [1]
Advanced—500µl, wash ratio=2 (as default).
3. Mix for 150.00 minutes at 450rpm
Note: the time can be changed. 90-150minutes preferred.
4. Empty CleavageBlock for 3.00 minutes
- 5.
6. Transfer 4000µl from CocktailSolution[1][] to CleavageBlock[1-2] using DCM [1]
Advanced—500µl, wash ratio=2.
7. Wait for 5.00 minutes
8. Empty CleavageBlock for 3.00 minutes
- 9.

Note: there should no mix step in the file after second flush, otherwise will cause contamination of the products.

FLU-2,CHM:

1. Flush Arm1 with DMF[1] and DMF[2]
- 2.
3. Repeat from step 1, 2 times
- 4.

Question: Each chemfile, there is EDIT→[SCALE]→ENTER SCALE FACTOR NUMBER? The scale factor will scale all volumes entered in the chemfile. For example if you wrote the chemfile with a volume of 1ml for a dispense, you could enter a scale factor of 4 and it would change the volume to 4ml.

[illegible]

[illegible]

Racks needed to be calibrated: RVBlock, Monomer, Piperidine, DIC, CleavageBlock, Cocktail Solution.

Note: all chemicals should have the same properties (e.g. MW) as what we have, otherwise, need to recalibrate.

Calibrate Arm 1:

Coordinates	X	Y	Z
Reference Pos	720	492	688
Cleaner Station1	20	102	2340
Waste Station1	20	170	2340

<i>RVBlock</i>	Diameter=20mm		X	Y
No. of Col.	5	Arm 1		
No. of Rows	8	Straight		
Vessel cross section area	314.16 mm ²	First Corner Pos	2344	2036
		Second Corner Pos	1831	634
		Z-max Pos.	2820	
		Z-Start Pos	2615	
		Z-Dispense Pos	2615	
		Z-Travel Pos	1	
		Aspiration Speed	18	
		Dispense Speed	15	
Define Vessels: activate used vessels and deactivate unused vessels.		Advanced: Can Empty. Recommend to use valve Front Chamber (3), deactivate unused valve chambers.		

Z-Max Pos: calibrate the probe 2-3 mm below the septa.

Z-Start Pos: the probe must **below the innermost septum** but above the level of the liquid.

Z-Dispense Pos: is the distance the probe will descend when delivering liquid. The probe must be 2-3 mm below the innermost septum.

Z-Travel Pos: is the position of the probe as the arm travels point to point on the tabletop. This should always be left to its default setting.

<u>Monomer</u>	Diameter=27mm		X	Y
No. of Col.	4	Arm 1		
No. of Rows	9	Straight		
Vessel cross section area	572.56 mm²	First Corner Pos	934	1985
		Second Corner Pos	493	128
		Z-max Pos.	2760	
		Z-Start Pos	1065	
		Z-Dispense Pos	1065	
		Z-Travel Pos	1	
		Aspiration Speed	18	
		Dispense Speed	17	
**Define Vessels (Important!!) Activate used vessels and deactivate unused vessels. Let the vessel 36 active.		Advanced: Not required.		

****IMPORTANT: Each time, you need to check and define your AAs. If not in the list, you need to calculate the MW. Different protecting groups will change the result.**

Z-Max Pos: calibrate the probe 2-3 mm from the bottom. **Note:** Z-Max high enough that the probe will always draw from the cylindrical region of the bottle.

Z-Start Pos: the probe must **below the innermost septum but above the level of the liquid** (Note: the bottle can not be too full).

<u>CleavageBlock</u>			X	Y
No. of Col.	5	Arm 1		
No. of Rows	8	Straight		
Vessel cross section area	314.16 mm²	First Corner Pos	2344	2036
		Second Corner Pos	1831	634
		Z-max Pos.	1826	

		Z-Start Pos	1768
		Z-Dispense Pos	1768
		Z-Travel Pos	1
		Aspiration Speed	20
		Dispense Speed	17
Define Vessels: activate used vessels and deactivate unused vessels.		Advanced: Can Empty. Recommend to use valve Front Chamber (3), deactivate unused valve chambers.	

<u>Piperidine</u>	Diameter=90 mm		X	Y
No. of Col.	1	Arm 1		
No. of Rows	1	Straight		
Vessel cross section area	6361.73 mm²	First Corner Position	1312	1998
		Second Corner Pos	1312	1998
		Z-max Pos.	2928	
		Z-Start Pos	1629	
		Z-Dispense Pos	1751	
		Z-Travel Pos	1	
		Aspiration Speed	18	
		Dispense Speed	17	
*Define Vessels (Important!!)		Advanced: Not required.		

*Define:

Chemistry Database	Reagent 1
Vessel No.	1
Full Name	Piperidine
Abbreviation	PIP
MW	85.20 g/mol
Density	1.45 g/ml

<u>DIC</u>	Diameter=90mm		X	Y
No. of Col.	1	Arm 1		
No. of Rows	1	Straight		
Vessel cross section area	6361.73 mm ²	First Corner Position	70	1014
		Second Corner Pos	70	1014
		Z-max Pos.	3023	
		Z-Start Pos	1554	
		Z-Dispense Pos	1719	
		Z-Travel Pos	1	
		Aspiration Speed	18	
		Dispense Speed	17	
**Define Vessels (Important!!)		Advanced: Not required.		

**Define:

Chemistry Database	Reagent 1
Vessel No.	1
Full Name	1,3-Diiso
Abbreviation	DIC
MW	126.20 g/mol
Density	0.81 g/ml

<u>Cocktail Solution</u>	Diameter=27mm		X	Y
No. of Col.	1	Arm 1		
No. of Rows	1	Straight		
Vessel cross section area	572.56 mm²	First Corner Pos	493	128
		Second Corner Pos	493	128
		Z-max Pos.	2760	
		Z-Start Pos	1065	

PART III: Synthesis Preparation:

Calculation: Here just show an example, yours may be different.

H-Val(V)-Gln(Q)-Ala(A)-Ala(A)-Ile(I)-Asp(D)- Ile(I)-Gly(G)-Tyr(Y)- OH

Rack	Vessel	Reagent	MW	Dens	M	Vol, ml	Dvol, ml	Total ml	Vol,	Grams
		DMF(1)				774.8		AA Total Vol: <u>46.4</u>		
		MEOH(2)				176.0				
AA	8	Gly	297.3			5.04	0.76	5.8		0.86
	10	Ile	353.4			10.08	1.52	11.6		2.04
	4	Asp	411.5			5.04	0.76	5.8		1.2
(AlaH2O)	1	AlaH2O	329.3			10.08	1.52	11.6		1.9
	7	Gln	610.7			5.04	0.76	5.8		1.78
	20	Val	339.4			5.04	0.76	5.8		0.98
PIP	1	PIP	85.2			68	30	98		4.17
DIC	1	DIC	126.2			24.16	30	54.16		3.42

(20-31: dep-2.cht, wash-1, couple, wash-2)

1. Print the calibration result, and check whether we have enough required AAs and reagents.
Note: if the protecting group is different from the calibration or the molecular formula is different, some modifications need to be considered.
2. Calculate the total volume of 0.5 M HOBt needed by adding Vol and Dvol columns together from all AAs listed on the calculation sheet.
Here: 46.4 ml
3. Prepare the **0.5 M HOBt solution in DMF (volume from above 2)**. For example, 38.27gHOBt+DMF⇒500ml (0.5M HOBt)}

HOBT.H ₂ O, g	Total volume, ml
7.65	100
15.31	200
22.96	300
30.62	400
38.27	500
45.92	600

7

4. Prepare a **25vol.% Piperidine in DMF** solution by adding the Vol and Dvol listed in calculations as the total volume of the solution, and place Piperidine in the calibrated location on the table top (**right front**). For example, 200ml Piperidine+600ml DMF⇒800ml (>Vol+DVol) 25vol% Piperidine in DMF (for deprotection).

Here: Total volume—98 ml, so prepare 160 ml (PIP: 40 + DMF: 120ml = Total 160ml)

Reagent Rack:

	36	35	34	33	
	Reagent K	...	...	...	
	32	31	30	29	
	28	27	26	25	
	24	23	22	21	
	20	19	18	17	
0.5 M DIC in DMF	Val (V)	Tyr (Y)	Trp (W)	Thr (T)	Reagent 1 (TFA)
DIC	16	15	14	13	
	Ser (S)	Pro (P)	Phe (F)	Met (M)	HOBt
	12	11	10	9	
	Lys (K)	Leu (L)	Ile (I)	His (H)	
	8	7	6	5	
	Gly (G)	Gln (Q)	Glu (E)	Cys (C)	PIP 25% PIP in DMF
	4	3	2	1	
	Asp (D)	Asn (N)	Arg (R)	Ala (A)	

5. Prepare **0.5 M DIC in DMF** solution by adding the Vol and Dvol listed in the calculations as the total volume of the solution, and place DIC in the calibrated location on the table top (**left behind**). For example, 23.23ml DIC+DMF⇒300ml (>Vol+DVol) 0.5 M DIC in DMF. For example, Dilute 7.74 ml DIC with DMF to total 100 ml.

Here: Total Volume—54.16 ml, so prepare 100 ml.

DMF	DIC	
100 ml	7.74 ml	6.31 g
200 ml	15.48 ml	12.62 g
300 ml	23.23 ml	18.93 g

Check: place the **monomer bottles** in the corresponding location in the monomer rack and place the top plate. Place **DIC (left behind)** and **Piperidine (right front)** in the calibrated locations on the tabletop (Extra DIC and PIP solutions store in hood, extra AA solutions store in -20°C fridge).

6. Check the *waste Carboy bottles* (will they overflow?), Fill **DMF for Dilutor 1** and **Dilutor 2** (may need to re-fill during synthesis). Important: should have enough DMF at hand (see calculation table for DMF amount).

7. Weigh all AAs into labeled monomer bottles.

Rack	Vessel	Reagent	MW	Dens	M	Vol, ml	Dvol, ml	Total ml	Vol,	Grams
		DMF(1)				774.8		AA Total Vol: 46.4		
		MEOH(2)				176.0				
AA	8	Gly	297.3			5.04	0.76		5.8	0.86
	10	Ile	353.4			10.08	1.52		11.6	2.04
	4	Asp	411.5			5.04	0.76		5.8	1.2
(AlaH ₂ O)	1	AlaH ₂ O	329.3			10.08	1.52		11.6	1.9
	7	Gln	610.7			5.04	0.76		5.8	1.78
	20	Val	339.4			5.04	0.76		5.8	0.98
PIP	1	PIP	85.2			68	30		98	4.17
DIC	1	DIC	126.2			24.16	30		54.16	3.42

- Reactors: 40 reactor positions

40	39	38	37	36
...	...	...	...	...
5	4	3	2	1
...				

- RUN:**

- Note: In case of an unexpected interruption of Synthesis, first write down the **Synthesis Position** (or go to the log, refer to Manual p100).

After synthesis, recommend to go ahead with PART IV as soon as possible. The cleavage/deprotection should be done as quickly as possible to minimize the exposure of peptide to the cleavage reagent. Therefore, *I guess*, after finish, cleavage and washing should be done continuously.

[illegible]

PART IV: Preparation for Cleavage

- Note: this program will not automatically start after the Synthesis, need to initiate.**

- Note: Ensure the cocktail solution is in the AA rack, monomer bottle position 36.

[illegible]

Preventive maintenance: at the conclusion of each synthesis

- Routine reactor cleaning procedure:**

- Note:

- Plunger overload:** recommended

- ### Plunger overload (Procedure: Kranthi7):

- [illegible]