



**Certificate of Equipment Training
Steris Amsco Lab 250
LV 250 Laboratory Steam Sterilizer
Department of Biochemistry, Room 3130J**

Name: _____

Job title: _____

Department: _____

Principal investigator: _____

Email: _____ **Phone:** _____

- I understand that high-temperature/high-pressure steam is used by the sterilizer to operate.
- I understand that my safety depends on the correct use of personal protective equipment: heat resistant gloves, goggles, lab coats and safety compliant footwear.
- I understand that sterilization occurs at temperatures between 120 and 130 C and that any liquid or plastic-ware must be contained in either a metal pan or a heat resistant plastic tray during the sterilization process.
- I have been instructed on the safe and correct use of the Steris Amsco Lab 250 sterilizer.
- I will fill out the logbook completely every time I use this instrument.
- I will keep the access door to room 3130J locked at all times
- I will not use this equipment if it is malfunctioning and I will report any malfunction to Paolo Fagone (Rm 3120, Tel. 304-293-7728; e-mail: pafagone@hsc.wvu.edu).

Trainee Signature/Date: _____

Principal investigator/Date: _____

Trainer Signature/Date: _____